



# Anesthesia apparatus Anes II for Mice & Rats



Easy to use and inexpensive type has just been released. It is designed for mice and rats. Omitted unnecessary options or functions and reduced the weight and cost. Just put the Isoflurane and adjust the concentration by the dial.

**Low Price! Please inquire**

Contact me :cbm03765@nifty.com

## Specification

Model	Anes II YS-18
Concentration	0~4.5%
Airflow	250ml/min accuracy $\pm 10\%$ *500ml/min type is also available.
Volume	Approx. 15ml
Temperature to use	20C-26C
Size, Weight	W180mm D95mm H100mm, 980g
Power	AC100V 0.5W
Standard accessories	Silicone mask for mice/rats (1ea), Tubing 5 x 7mm 1meter, 16G needle for filling Isoflurane, Extra gas absorbing filter x 1

**BioMachinery Co., Ltd.**

1-9-1-111 Sakae, Funabashi Chiba. 273-0018, Japan Tel (81) 47-432-0091 Fax (81) 47-431-7979